JEILMEDICAL

Dual Top[®] ANCHOR SYSTEM

Instruction For Use

- Dual Top Anchor System
- Bone Plate/Bone Screw/Sterile Bone Screw

Issue Date: 2024. 04

Print Date:

#497-41002-04

Table of Content

English	03
Deutsch	14
Français	26
Português	38
Polski	50
Ελληνικά	62
Nederlands	74
Čeština	86
Română	98
Magyar	110
Slovenčina	122
Svenska	134
Norsk	146
Italiano	158
Türkçe	170
Español	182

[INTENDED PURPOSE]

Dual Top Anchor System is used to provide a fixed anchorage point for attachment of orthodontic appliances to facilitate the orthodontic movement of teeth. It is used temporarily and is removed after orthodontic treatment has been completed. For use on patients aged 12 or above who need orthodontic treatment.

[MATERIAL]

Bone Screw : Titanium Alloy (ASTM F 136)

Bone Plate : Titanium (ASTM F 67)

[INTENDED USERS]

Dual Top Anchor System bone screw and bone plate are intended to be used by dentists or clinicians who are trained, experienced, and familiar with the use of orthodontic anchorage techniques. The dentist must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.

[PATIENT TARGET GROUP]

The device is to be used with consider intended use, indications, contraindications and patient's bone anatomy and health condition. The intended patient should be aged 12 or above.

- Indication : Malocclusion Class I, II, III

[CLINICAL BENEFIT]

- Stabilization of orthodontic anchorage and relief of malocclusion symptoms are expected to be achieved from successful treatment in patients with malocclusion during orthodontic treatment.
- SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.
 - Eudamed public website (URL): <https://ec.europa.eu/tools/eudamed>
 - Basic UDI DI: 08800089410033L2
 - UDI Carrier

	UDI Information	AI Information
UDI-DI	(01)0880XXXXXXXXXX	(01) GTIN
UDI-PI	Non-Sterile product (11)YYMMDD(10)JYYMMDD-LXXX	(11) Production date (10) Lot No.
	Sterile Product (17)YYMMDD(10)JYYMMDD-LXXX	(17) Expiration date (10) Lot No.

[PERFORMANCE CHARACTERISTICS OF THE DEVICE]

- The performance and safety of the Dual Top Anchor System bone screw and bone plate have been verified and that represent state of the art medical devices for orthodontic anchorage techniques on to mandible and maxilla.

[DESCRIPTION]

The Dual Top Anchor System was designed to provide a firm anchorage in intraorally mandible and maxilla for an orthodontic treatment. Bone Plates are manufactured from titanium (ASTM F67) and Bone Screw are manufactured from titanium alloy (ASTM F136). The Dual Top Anchor System is designed to apply to various orthodontic tools and diverse enough to satisfy various clinical cases.

[DIRECTIONS FOR USE]

After providing the patients with a thorough explanation on the necessity, benefits and shortcomings of using the Dual Top Anchor System, dentist must obtain the patient's consent prior to use.

1. BEFORE IMPLANTATION

- It is necessary to check if the patient has any allergies.
- Determining the implant location

: Choose the implantation site according to the treatment objective and the quality and quantity of bone. Use any metal device and x-ray image to confirm that the implantation site allows adequate clearance from the tooth roots and nerves.

3) Anesthesia in the designated area for screw insertion

: Disinfect the implant site with an antiseptic solution then inject the anesthetic. Make sure that the anesthesia is limited to the peripheral soft tissue.

2. AFTER IMPLANTATION

- 1) Proceed with the subsequent orthodontic procedure according to the dentist's judgment.
- 2) The recovery period can be extended if adequate fixation force is not achieved due to poor bone quality.
- 3) During the treatment period, at patient follow-up visits, orthodontists should check for hygiene, possibility of latent inflammation, and stability of the screws.
- 4) If necessary, provide the following information to the patient.
 - ① Avoid hard or tough food for about a week after the procedure and refrain from touching it with your hands or tongue
 - ② Avoid drinking and smoking for 1 to 2 weeks after the procedure
 - ③ To prevent inflammation, use a soft toothbrush to keep the area around the screw clean
 - ④ If you are shocked around your face or mouth while using the device, or if the device is shaken, visit the hospital.

3. AFTER DEVICE REMOVAL

- 1) If necessary, provide the following information to the patient.

① Maintain oral hygiene to prevent inflammation.

4. IMPLANTATION PROCEDURES

1) Bone Screw

- ① Select the screw based on the implantation location and the condition of the patient's bone quality.
- ② Attach the selected screw to the compatible driver. The driver should be perpendicular to the screw.
- ③ Place the screw vertically on the surface of the mucosa, and then turn the handle clockwise while pressing it so that the tip of the screw can dig a hole in the cortical bone. Turn the driver slowly and with caution in order to prevent the screw axis from shaking.
- ④ Stop inserting when the head at the lower end of the screw meets the mucosa, and then remove the driver from the screw.

2) Bone Plate

- ① After the diagnosis was concluded about the orthodontic treatment, the dentist selects the surgical site. However, if an incision is required at the location where the plate will be placed, it will be determined by the dentist's judgment.
- ② Select the screw and plate based on the implantation location and the condition of the patient's bone quality. Plate may need to be reshaped and/or trimmed to conform to the patient's bony anatomy.
- ③ Attach the selected screw to the compatible driver(20-AT Series). The driver should be perpendicular to the screw.
- ④ While holding the plate in the desired location, and place the screw vertically on the surface then rotate clockwise while pressing it to fixate the plate and screw in the surgical site.

5. BONE PLATE & BONE SCREW REMOVAL

It is recommended that the implant be removed after use for up to 24 months, which is the period of completion of orthodontic treatment, unless medical attention is required. However, the removal period may be extended if there is a special opinion of the dentist. Remove the orthodontic appliances associated with bone plate/bone screw and use a screwdriver to turn the screw counterclockwise to remove the screw. Remove all the screws, remove the plate and finish.

6. RECOMMEND USAGE SITE

• It can be used in any part of the oral cavity, but it is recommended that it is selected as a fixed anchorage point for orthodontic treatment by selectively placed in the maxillary and mandible parts under the decision of the doctor.

- 1) Maxilla
- 2) Mandible

7. ORTHODONTIC APPLIANCE

※ Orthodontic appliance is available when they are compatible with the product regardless of brand.

- 1) Bone Screw

Type	Orthodontic Appliance	Driver shaft
JA Series	Rubber band	MJ, MD Series
	Arch wire	
JB Series	Rubber band	JB Series
	NiT Spring	
JD Series	Arch wire	MJ Series
G1 Series	Arch wire	GD Series
G2 Series	Rubber band	MJ, GD Series
	Arch wire	
JF Series	Rubber band	MJ Series
	Arch wire	
JS Series	Rubber band	MJ, MD Series
	Arch wire	

(Size : Rubber band – 1/8"~3/8", Arch wire – 0.012~0.018", 0.018X0.025"~0.021X0.025", Ni-Ti Spring : Ø1.8 ~ Ø2.0mm)

2) Bone Plate

Type	Orthodontic Appliance	Screw
AP-T series	Rubber band Arch wire	20-AT Series
AP-L series		
AP-S series		

(Size : Rubber band – 1/8"~3/8", Arch wire – 0.018X0.025"~0.021X0.025") Only designated screw is available for Plate.

[WARNING]

1. Selection of plate and screw size must be carefully considered by the operating dentist and should take into consideration the quality of bone, bone type, functional loads exerted on bone(s), and post-operative patient compliance.
2. Avoid bending of the screw hole, sharp bending and multiple bending as it increases the risk of implant breakage. bending as it increases the risk of implant breakage.
3. The dentist must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.
4. Take care to use of an orthodontic appliance that can cause excessive, abnormal stress on all devices.
5. Patients whose device have been implanted for more than 24 months need to continuous observation of complications.
6. Use of an undersized plate or screw in areas of high functional stress may lead to implant fracture and failure.

[CAUTION]

1. Fixing the product incompletely may cause damage or breakage of the product, and the patient may feel a foreign body sensation.
2. It is at the physician's discretion to evaluate the patient's medical conditions and select a fixation device most appropriate for the individual patient.
3. Responsibility for proper selection of patients, adequate training, experience in the choice and placement of plate & screw and the decision to leave or remove plate and screw postoperatively, rests with the dentist.
4. When engaging the screw, axial pressure of the driver shaft into the screw head must be adequately applied to ensure that the blade is fully inserted into the screw head. This results in proper axial alignment and full contact between driver shaft and screw, minimizing the risk of round-out. Otherwise, there will be an increased risk of mechanical damage to the implant or the driver shaft.
5. Inserting the screw with excessive force may cause damage or breakage of the product, especially if the bone density of the patient's is too high or the bone thickness is too thick.
6. If contouring is necessary, the surgeon should avoid bending the device at a screw hole.
7. Be aware that the patient may be injured if the device is swallowed.

[PRECAUTIONS]

1. Before use, check that the packaging & product are not damaged.
2. The dentist should inform the patient of risks associated with surgery.
3. Before use, check that the specifications of the product matches those on the label.
4. Check the brightness of the operating room is appropriate.
5. User should read the instructions before use. It is the professional's responsibility to use the Jeilmedical products according to the instructions for use.
6. Do not use products that fall on the floor. This product should be sterilized before use.
7. Check the compatibility of instruments before use.

[CONTRAINDICATION]

1. Not for use in cases of active or suspected infection or in patients previously sensitized to titanium.
2. Not for use in patients exhibiting disorders which would cause the patient to ignore the limitations of rigid fixation plate and screw implants.
3. Not for use in patients with children (<12 years), uncontrollable diabetes, excessive smokers or alcoholics.
4. Do not implanted in the tooth roots, vessels, nerves and sinuses.
5. Not for use in patients with poor density or poor bones in implantation site.
6. Not for use in patients with pregnant/breastfeeding women.

[SINGLE-USE ONLY]

The plates and screws are disposable. Do not reuse. Any implant should never be reprocessed or reused after contact with blood, body tissues or fluids, but rather should be disposed. Reprocess or reuse of device may potentially result in serious patient harm such as infection. Jeil Medical Corporation does not take any responsibility for all problems caused by reuse.

[SIDE EFFECTS AND RESIDUAL RISK]

Residual risks, side effects and adverse events may occur in all surgical procedures. Dentists should operate carefully to prevent problems and take appropriate action for the patient after surgery. Undesirable side effects(allergy) may also occur.

Many reactions can occur and the most common are as follows:

- Infection
- Inflammation(Soft tissue inflammation, Gingival inflammation, Ulceration, Gingivitis or buccal tissue irritation, Hypertrophy in the gingival tissues)
- Root resorption
- Tooth damage
- Perforation(Perforation of the maxillary sinus and nasal cavity walls.)
- Bleeding(Reinstallation, Implant loss, fracture)
- Swelling
- Discomfort(Implant mobility and loosening)
- Pain(Distress)
- Growth disturbances
- Allergy(Hypersensitivity reactions, Thermal sensitivity, Skin irritation, etc.)

[MR SAFETY INFORMATION]

MR Conditional

1. Radio Frequency induced heating according to ASTM F 2182-19e2

Radio frequency induced heating evaluation of a worst case device lead to temperature rises of 5.45°C(1st probe), 0.51°C(2nd probe), 10.53°C(3rd probe) at 1.5 T and 2.82°C(1st probe), 1.00°C(2nd probe), 3.44°C(3rd probe) at 3.0 T under MRI Conditions using FDA approved MR Coils.

2. Torque and Deflection angle according to ASTM F2213-17, ASTM F2052-15

Deflection angle test was repeated three(3) times and it was found that the samples were moved 0° at the highest condition of MRI coil. The torque testing was repeated three(3) times at the isocenter and the samples moved 0 degree in all direction. Therefore quantitative scale was 0. No torque was detected under 1.5T MRI.

3. Image artifacts according to ASTM F 2119-07

Under T1 SE and GRE pulse sequence, product extended maximum 4.64mm artifact size compared to its original shape.

The above mentioned test relies on nonclinical testing. The actual temperature rise in the patient will depend on a variety of factors beyond the SAR and time of RF application. The patient using the device can be scanned safely under following conditions:

- Static magnetic field of 1.5T MRI, 3.0T MRI, only
- Maximum spatial gradient magnetic field of 720 Guass/cm(extrapolated) or less
- Maximum MR system reported, whole body averaged SAR of 2 W/kg

[STERIZATION]

• For Non-Sterile

Dual Top Anchor System must be sterilized before use. Following instructions are recommended;

- ① Place the devices in the appropriate block (kit or tray) using forceps and/or powder-free gloves to avoid contamination and any other negative effect on the surface of device.
- ② Wrap the block with a surgical drape.
- ③ Sterilize in the autoclave validated and maintained in accordance with ISO 17665 and ANSI AAMI ST79. Following parameters are validated in accordance with ISO 17665-1 and recommended for sterilization;

Type	Temperature	Exposure Time	Load Characteristics	Drying Time
Gravity	132°C	15 min.	Wrapped*	30 min.
Pre-vacuum	132°C	4 min.	Wrapped*	30 min.

* In the case of load characteristics, we recommend usage of an FDA cleared wrap to ensure that the device is actually sterile prior to implantation.

• For Sterile

Dual Top Anchor System Sterile Bone Screw is supplied by gamma sterilization.

Use-by date : Refer to label

[STORAGE AND HANDLING CONDITION]

Careful handling and storage of the product is required. Store in a dry place at 10~22°C temperature and 10~30% humidity.

[DISPOSAL]

Any implant that has been contaminated by blood, body tissues or fluids, should never be used again and should be handled according to hospital protocol. Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

[COMBINATION DEVICE]

Only compatibility between devices provided by Jeil medical coporation and Orthodontic Appliance of the size specified on the IFU has been confirmed (Driver shaft, Hand body, Rubber band, Arch wire, NiTi spring). Compatibility with other manufacturer's devices have not been tested and there is no liability in such instances.

[INCIDENT NOTIFICATION]

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

[PACKAGING]

Dual Top Anchor System are packaged individually or set. This device manufactured and sold by Jeil Medical Corporation.

[SYMBOL DESCRIPTIONS]

<i>Symbol</i>	<i>Description</i>	<i>Symbol</i>	<i>Description</i>
	Catalogue number		Batch code
	Do not reuse		Do not resterilize
	Date of manufacture		Manufacturer
	Non-Sterile		Use-by Date
	Sterilized using irradiation		Caution
	Consult instruction for use		Do not use if package is damaged
	Authorized representative in the European Community		EMERGO EUROPE Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands
	Unique device identifier		Medical device
	Temperature limit		Humidity limitation
	Single sterile barrier system with protective packaging inside		

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