INDICATIONS :

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INDICATIONS :

The Herbst-Amoric Orthosis (HAA) is a component of the removable mandibular advancement appliance.

The appliance forces the mandible to move forward, enabling correction of tooth position and increasing oropharyngeal and/or hypopharyngeal air space.

The OHA appliance is a single-use medical device for a single patient.

It is indicated for the correction of malocclusion (Class II), the treatment or alleviation of mild to severe Obstructive Sleep Apnea/Hypopnea Syndrome, and the treatment of certain temporomandibular joint pathologies.

The OHA unit is intended for use by dental technicians from specialized laboratories. The device is made available to end-users (qualified dentists/orthodontists) who have the necessary skills for safe handling, and who are responsible for patient treatment and follow-up.

Remarks :

Adjustments for patient's comfort are made by shortening the length of the tube or by placing crimped mini-shims on the rod. It is also possible to add vertical elastic rods on the gutters to force the patient to close the mouth during their sleep.

PRECAUTIONS FOR USE :

The OHA device is supplied NON-STERILE. Cleaning and sterilization in compliance with local regulations are recommended before each first use to avoid any risk of cross-contamination.

Healthcare professionals are responsible for providing patients with the information they need to use the device, taking into account the following recommendations:

- Always check the unit for cleanliness, abnormal signs or damage before placing it in the mouth.
- Routinely rinse the appliance before placing it in the mouth and after removing it from the mouth.
- Respect the recommended wearing time.
- Use oral hygiene products without fluoride components (risk of corrosion).

Storage conditions : Store at room temperature in a dry, well-ventilated place away from direct sunlight.

Accessories: lingual shields for thermoformed trays can be used during treatment.

CONTRAINDICATIONS :

- Do not use ethanol for cleaning
- Do not soak device in hot water
- Do not use these braces in any of the following situations:
 - Allergy to NICKEL
 - Missing teeth
 - Presence of damaged or decayed teeth, periodontal disease with dislocation, arthritis and ankylosis of the temporomaxillary joints.

Healthcare professionals must inform patients of any possible adverse effects.

Warning: certain side effects may occur during treatment:

- normocclusion, overbite or class III malocclusion may be exacerbated
- pain in the temporomandibular joints and masticatory muscles
- dental or jaw joint discomfort.
- dry mouth, abnormal occlusion, occlusal problems, tooth mobility, excessive salivation.

IMPORTANT

No liability can be accepted for products which have been modified or used for purposes other than those for which they were intended, or which have been used inappropriately. All the instructions specified above have been validated as suitable for the treatment of a medical device for the purpose of reuse.

Any adverse reaction linked to medical devices should be notified to the manufacturer, and any serious incident should be reported without delay to the competent authority of the Member State in which the patients and practitioners are based.

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