INDICATIONS

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INDICATIONS

The orthodontic instruments at your disposal allow you to cut, hold wires or small parts, and form bands, ligatures, etc. These are medical devices that come into contact with mucous membranes and skin. They are intended to be used by health professionals who must have sufficient knowledge to guarantee their proper handling and sterility at each use.

PRECAUTIONS FOR USE AND STORAGE

The legal regulations in force in your country as well as the hygiene regulations of the doctor's surgery or hospital must be observed.

Improper handling and maintenance, as well as improper use, can lead to premature wear of the orthodontic instruments in your possession.

The user is responsible for the sterility of the instruments and must ensure that validated cleaning, disinfection and sterilization processes are used.

- **Before each 1^{ère} use:** New orthodontic instruments must be thoroughly cleaned, disinfected and sterilized.
- **As a routine before each use:**
 - . Examine each instrument thoroughly to ensure that it is not broken, cracked or defective, particularly at the level of the tips of the jaws, and pins.
 - . Check the condition of the locking devices and moving parts.
 - . Check that the size limit of the wire to be cut is respected, if applicable.
- **Protection of personnel and third parties:** Contaminated instruments should be disinfected as soon as possible during a treatment process to ensure maximum safety for the personnel handling them. When handling all used and contaminated instruments, it is necessary to wear protective gloves that comply with the requirements of Regulation (EU) 2016 / 425 on Personal Protective Equipment (PPE).
- **Cleaning, disinfection and sterilization processes:** Only validated processes are used. Effective cleaning and disinfection are prerequisites for the correct sterilization of orthodontic instruments. In addition, the sterilization equipment must be regularly maintained and checked. The validated parameters for the cleaning, disinfection and sterilization cycles must also be monitored regularly.
- **Storage condition:** In rooms with low humidity to avoid condensation.

CONTRAINdications

- . Do not use damaged orthodontic instruments.
- . Do not repair yourself.
- . Do not have untrained personnel service or repair orthodontic instruments.
- . Do not use NaCl solution or wire brushes or steel wool and/or abrasive drying cloths/papers.
- . Do not use lubricants with additives

INSTRUMENT PROCESSING

- **Pretreatment :**
 - Remove impurities immediately after use, preferably within 2 hours
 - Use only an approved solution of a cleaning agent/disinfectant mixture (with CE mark) without protein binding effect.

Important: In case of mixing, it is imperative to respect the recommendations of the manufacturers of the products.

Disinfectants used during pretreatment are for personal safety only and are not a substitute for subsequent instrument disinfection.

- Do not place orthodontic instruments in NaCl solutions, (risk of corrosion)
 - Avoid overloading instrument baskets and wash trays
 - Always process jointed instruments in the open position.
 - Use a soft brush to manually remove coarse impurities, not wire brushes or steel wool (risk of scratching)
-
- **Cleaning:** pH neutral or alkaline cleaning products
 - **Automatic cleaning :**
 - 1 - Place the jointed instruments in the device so that the joints remain open and water cannot stagnate
 - 2 - Wash and rinse according to the manufacturer's instructions.
 - 3 - Remove the instruments at the end of the cleaning cycle and visually check the cleanliness of the jaws, shafts and handles; if necessary, repeat the automatic cycle or perform a manual cleaning.
 - **Manual cleaning :**
 - 1 - Thoroughly clean surface dirt from the instrument
 - 2 - Spread the cleaning solution with a brush on all surfaces
 - 3 - Clean the instruments with joint in both open and closed position.
 - 4 - Rinse the orthodontic instrument under running water. The water should flow through the cannulas, and the blind holes should be filled and emptied several times.
 - **Ultrasonic tank :**
Its use is recommended in the cleaning process; the immersion time should not exceed the time recommended by the manufacturer of the detergent used.
 - **Disinfection :** This follows the cleaning process; in the case of automatic cleaning, a 10-minute thermal disinfection at 93°C can be performed according to the equipment manufacturer's instructions. Disinfection solutions must be used according to the manufacturer's instructions.
 - **Final rinse:** Use demineralized water whenever possible. This will prevent stains, deposits and corrosion of the instrument.
 - **Drying:** In manual mode, using non-abrasive paper/cloth to avoid scratching.
If drying is an integral part of the cleaning/disinfection cycle, a temperature of 93°C should not be exceeded.
 - **Sterilization:**
The recommended sterilization method is wet heat sterilization according to ISO 17665-1: 2006
Sterilizers and sterilization processes must comply with current standards, and more precisely to the requirements described in the standard EN NF 13060 +A1 - Small steam sterilizers
 - Make sure to respect the maximum load allowed by the equipment manufacturer during a sterilization cycle (case of instrument pooling)
 - Make sure to respect the parameters of sterilization time and temperature defined as follows:
15 min at T=121°C or 3 min at T=134°C

CAUTION: rapid sterilization leads to faster wear of the instrument

INSTRUMENT MAINTENANCE

- **Before sterilization :**
 - Lubricate with a small amount of silicone spray or high quality medical white oil on joints and moving parts.
 - Visually check for cracks, holes, microporosities or damage and remove dull or damaged instruments if necessary
 - Check the mobility of the jointed instruments (avoid too much play) and the correct operation of the locking mechanisms (locking wheel)

• Bagging/packaging before sterilization :

- Use bagging equipment that complies with standards; the bag must be large enough not to be stretched.
- Place orthodontic instruments in the trays provided or place them on universal sterilization trays.
- Protect the jaws if necessary.
- Use an appropriate process to pack the trays.

• Recommended cleaning products :

- ECI0150 : White oil medical instrument 50 ml
- ECM1070 : IP-Spray Instruments 400 ml

• Warranty:

- The warranty is excluded in the case of products that have been modified or used outside of their intended use or in an inappropriate manner
- The manufacturer's warranty no longer applies in the event of repair or modification by anyone other than the manufacturer.
- All of the instructions specified above have been validated as appropriate for processing a medical device for reuse.

ATTENTION:

Returned instruments must be disinfected and sterilized, accompanied by a document attesting to the sterilization and the date it was performed.

IMPORTANT :

Any serious incident arising from the use of this device shall be reported to the manufacturer and to the regulatory authority of the Member State where the patient/user is established.

In compliance with the applicable statutory provisions and, above all, for the protection of our employees, we require a signed "Proof of decontamination" for each return. Please ensure that this "Proof of decontamination" is filled in and enclosed with each product return (complaint / repair / other reason for return) and package the product such that no risk of injury exists for our personnel in incoming goods inspection while unpacking.

Article no.: Article designation:

If possible provide information on

Delivery note N°: of [date] Batch N°:

Reason for return:.....

We hereby confirm that

- The medical device(s) enclosed with the letter has/have not come into contact with blood, tissue or other bodily substances/fluids and hygienic safety can be confirmed by signing.
- The medical device(s) enclosed with the letter has/have come into contact with blood, tissue or other bodily substances/fluids and have been disinfected, cleaned and sterilized according to the applicable hygiene requirements for medical devices and the manufacturer's specifications, which we confirm by signing.

Specifications on disinfection, cleaning and sterilization:

- Disinfection and cleaning have been performed manually
 - Disinfection and cleaning have been performed mechanically
 - Steam sterilization (3 min. at 134°C / 15 min. at 121 °C)
 - Other method (please specify)
- The medical device(s) enclosed with the letter could not be decontaminated!

Reason:.....

Stamp from the institution

Date

Signature (and full name)

IMPORTANT : Should we not receive this document or a comparable confirmation, we reserve the right to return the goods to the sender at their own expense.

SURVEY ON MAINTENANCE METHODS OF ORTHODONTIC INSTRUMENTS

The informations below are specific to procedures that were implementing in your office. They will also help us to give you an appropriate answer to the issue met. Thank you in advance for your cooperation,

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1 – CLEANING :	<input type="checkbox"/> YES	<input type="checkbox"/> NO	
<input type="checkbox"/> Manual			
<input type="checkbox"/> Automatic	Material : -----		
Time of cleaning: -----			
PRODUCT(S) : -----			
WATER :	<input type="checkbox"/> tap water	<input type="checkbox"/> demineralised water	<input type="checkbox"/> Other -----
RINSING :	<input type="checkbox"/> YES	<input type="checkbox"/> NO	
If yes, which kind of water ?: -----			
DRYING	YES <input type="checkbox"/>	NO <input type="checkbox"/>	
Tissue	<input type="checkbox"/>		
Compressed air	<input type="checkbox"/>		
others	<input type="checkbox"/> -----		

2 – DECONTAMINATION :	<input type="checkbox"/> YES	<input type="checkbox"/> NO	
<input type="checkbox"/> Manual			
<input type="checkbox"/> Automatic	Material : -----		
Time of decontamination : -----			
PRODUCT(S) : -----			
WATER :	<input type="checkbox"/> tap water	<input type="checkbox"/> demineralised water	<input type="checkbox"/> Other-----
RINSING :	<input type="checkbox"/> YES	<input type="checkbox"/> NO	
If yes, which kind of water ?: -----			
DRYING	YES <input type="checkbox"/>	NO <input type="checkbox"/>	
Tissue	<input type="checkbox"/>		
Compressed air	<input type="checkbox"/>		
others	<input type="checkbox"/> -----		

3 – STERILIZATION :	<input type="checkbox"/> YES	<input type="checkbox"/> NO
<input type="checkbox"/> Cold sterilization – product : -----		
<input type="checkbox"/> Autoclave	Time of cycle : ----- / T°maximum : -----	
Protection : <input type="checkbox"/> YES – which way (bags, oiling...) : -----		
<input type="checkbox"/> NO		

4 – LUBRIFICATION :	<input type="checkbox"/> YES	<input type="checkbox"/> NO
When ? ----- / Frequency ? -----		
Product(s) : -----		

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