

Distributed by



ortho·jel CE

ALGINATE POUR EMPREINTE

Directions For Use

For dental use only

1. Transfer the impression material from the bag to the container.
2. Invert the container two or three times to loosen any compacted material before using.
3. For best results mix with water at 23°C. Colder water retards setting, warmer water accelerates setting (see table below) :

Temp . (°C)	Mixing time (seconds)	Approx. Working time from start of mixing (seconds)	Approx. setting time from start of mixing (seconds)
21	55	95	125
23	50	90	120
25	45	85	110

4. Have patient rinse mouth out with warm water.
5. Lightly dip the large scoop into the powder, tap gently and remove any excess with a spatula.
6. for normal impressions, one large scoopful (16 g) is sufficient.
7. Add one measure of water (small scoop 39ml) and spatulate thoroughly until a creamy smooth mix is obtained. Ortho.Jel will be violet coloured during mixing.
8. Load the impression tray when the mix colour turns pink.
9. When the mix starts to turn white, insert the loaded tray into the patient's mouts.
10. Position the tray in the mouth, muscle trim and hold steady until set.
11. Remove the impression from the mouth carefully abd rinse it with cold water, ensuring that all saliva and debris is eradicated.
12. Where possible, alginate impressions should be cast immediately..
13. Alginate impressions should be stored and shipped only in a moist atmosphere of 100% (i.e. in a polyethylene bag along with a damp tissue) and cast within 48 hours of taking the impression. Gypsum products recommended for use with Ortho Jel are Moldana, Calestone, Fuji Rock and Velmix.
14. Do not store the set impression in water.
15. After each use, firmly close the lid of the container.
16. Store in a dry place at room temperature(23 ± 2°C / 76 ± 6°F)

Type 1 Fast Set Class A

Conforms to BS 4269: Part 2/ EN21563/ISO 1563/ADA 18/DIN 13909/AS1282

Manufactured in Scotland by

**WRIGHT HEALTH GROUP LTD,
Kingsway West, Dundee**