Prior to use, carefully read the instructions for use.

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GENERAL

In all cases

- Sandblast metal surfaces with alumina (aluminium oxide), and blow clean with air gun.

Note :

- Avoid contamination with oil, separator or moisture from fingers, etc.
- Dispense 2-3 drops of METALPRIMER II into a dappen dish or similar container. Using the brush provided, apply a thin layer to the bonding surface. Allow to dry.

Note :

- A second coat may be applied if necessary.

- Apply resin to metal as indicated below.

Note :

- Apply resin to metal immediately after the surface is prepared.

RECOMMENDED APPLICATIONS

- Resin cement to metal for metal inlays, crowns and bridges
 - Prepare metal surface as above and bond restoration with resin cement.
- Opaque resin to metal for composite resin veneers and crowns
 - Prepare surface as above and apply opaque (GC GRADIA) in usual manner. Cure as recommended by the manufacturer.
- Denture base to magnetic attachment
 - Prepare surface as above and bond attachment in usual manner using self curing acrylic resin (GC UNIFAST Trad).
- Denture base resin to metal base, clasp or attachment
 - Prepare bonding surface with alumina sandblast and blow clean with air gun.
 - Form a wax denture.
 - After flushing out the wax, clean the bonding surface using a steam cleaner. Remove wax completely using solvent.
 - Apply METALPRIMER II.
 - Pack and cure denture base in usual manner.
- In the mouth repair of damaged crown and bridge veneers
 - Place rubber dam.
 - Clean metal surface with carbide bur and prepare as above.
 - Special note :

Be sure to use rubber dam and avoid direct contact with oral tissue. In case of contact, immediately flush with water and seek medical attention.
- Repair of denture with metal base, clasp or attachment
 - Roughen metal bonding surface using carborundum point.
 - Prepare bonding surface as above.
 - Build up using self-cure resin in the usual manner.
- Bonding resin to metal orthodontic appliances
 - Prepare bonding surfaces as above.
 - Apply orthodontic resin to the appliance in the usual manner.

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C **E** 0086

GC.

MANUFACTURED by
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium
TEL: +32. 16. 39. 80. 50.

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900.

GC ASIA DENTAL PTE. LTD.
19 Loyalng Way, #06-27 Singapore 508724
TEL: +65 6546 7588

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EN Prior to use, carefully read the instructions for use.

METALPRIMER II

PRIMER FOR BONDING RESIN TO METAL

METALPRIMER II contains a special monomer (MEPS : thiophosphoric methacrylate) which promotes bonding of acrylic resins and composites to all types of metal.

RECOMMENDED USE

Bonding of acrylic and composite resins to all types of metal as set out below in the DIRECTIONS FOR USE.

NOTES AND CONTRAINDICATIONS

- For use only by dental professionals for the applications recommended.
- Contains methacrylate. Not to be used by operators, or on patients, with a history of allergy to methacrylates. In cases of allergic response, seek medical attention.
- In case of contact with skin or eyes, immediately flush with water and seek medical attention.
- Flammable. Avoid high temperatures, direct sunlight and sources of ignition. Use in a well ventilated area.
- Volatile. Seal container immediately after use.
- Assign the brush provided (No. 7) for use ONLY with METALPRIMER II.

DIRECTIONS FOR USE

- GENERAL**
In all cases
- Sandblast metal surfaces with alumina (aluminium oxide), and blow clean with air gun.**Note :**Avoid contamination with oil, separator or moisture from fingers, etc. 2.Dispense 2-3 drops of METALPRIMER II into a dappen dish or similar container. Using the brush provided, apply a thin layer to the bonding surface. Allow to dry.**Note :**A second coat may be applied if necessary. 3.Apply resin to metal as indicated below.**Note :**Apply resin to metal immediately after the surface is prepared.

RECOMMENDED APPLICATIONS

- Resin cement to metal for metal inlays, crowns and bridges
 - Prepare metal surface as above and bond restoration with resin cement.
- Opaque resin to metal for composite resin veneers and crowns
 - Prepare surface as above and apply opaque (GC GRADIA) in usual manner. Cure as recommended by the manufacturer.
- Denture base to magnetic attachment
 - Prepare surface as above and bond attachment in usual manner using self curing acrylic resin (GC UNIFAST Trad)
- Denture base resin to metal base, clasp or attachment
 - Prepare bonding surface with alumina sandblast and blow clean with air gun.
 - Form a wax denture.
 - After flushing out the wax, clean the bonding surface using a steam cleaner. Remove wax completely using solvent.
 - Apply METALPRIMER II.
 - Pack and cure denture base in usual manner.
- In the mouth repair of damaged crown and bridge veneers
 - Place rubber dam.
 - Clean metal surface with carbide bur and prepare as above.
 - Special note : In case of contact with rubber dam and avoid direct contact with oral tissue.
- In case of contact, immediately flush with water and seek medical attention.

- Repair of denture with metal base, clasp or attachment
 - Roughen metal bonding surface using carbondurum point.
 - Prepare bonding surface as above.
 - Build up using self-cure resin in the usual manner.
- Bonding resin to metal orthodontic appliances
 - Prepare bonding surfaces as above.
 - Apply orthodontic resin to the appliance in the usual manner.

PACKAGE

- bottle (5 mL) METALPRIMER II
- brush (No. 7)

STORAGE

Store in a cool place (4-25°C) out of direct sunlight. Shelf-life - 2 years from date of manufacture.

Last revised : 10/2008

MSK 0086

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Loven, Belgium TEL: +32 16 39 80 80.

PRINTED IN JAPAN

MANUFACTURED BY
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Loven, Belgium TEL: +32 16 39 80 80.

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RECOMMENDED APPLICATIONS

- Resin cement to metal for metal inlays, crowns and bridges
 - Prepare metal surface as above and bond restoration with resin cement.
- Opaque resin to metal for composite resin veneers and crowns
 - Prepare surface as above and apply opaque (GC GRADIA) in usual manner. Cure as recommended by the manufacturer.
- Denture base to magnetic attachment
 - Prepare surface as above and bond attachment in usual manner using self curing acrylic resin (GC UNIFAST Trad)
- Denture base resin to metal base, clasp or attachment
 - Prepare bonding surface with alumina sandblast and blow clean with air gun.
 - Form a wax denture.
 - After flushing out the wax, clean the bonding surface using a steam cleaner. Remove wax completely using solvent.
 - Apply METALPRIMER II.
 - Pack and cure denture base in usual manner.
- In the mouth repair of damaged crown and bridge veneers
 - Place rubber dam.
 - Clean metal surface with carbide bur and prepare as above.
 - Special note : In case of contact, immediately flush with water and seek medical attention.
- Repair of denture with metal base, clasp or attachment
 - Roughen metal bonding surface using carbondurum point.
 - Prepare bonding surface as above.
 - Build up using self-cure resin in the usual manner.
- Bonding resin to metal orthodontic appliances
 - Prepare bonding surfaces as above.
 - Apply orthodontic resin to the appliance in the usual manner.

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