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- **Fonction**

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
High pull head strap : ECJ3920/3920C

- **Function**

Class I medical device which, when connected to an orthodontic appliance, serves as a support for traction.
Maximum load capacity of 2.5 kg per side for the release force of the safety module.

- **Indications**

For use during orthodontic treatment, based on the doctor's recommendation.

- **Contraindications**


Allergic reactions when wearing the cap during treatment.

Material : 100% cotton


Wash at 40°C with a general detergent used for cotton.

Discoloration risk : never wash together with other color headcap.

The product must not be reused if it is damaged. Any serious incident involving the device should be reported to the manufacturer and to the competent authority of the country in which the user resides.

Usage

safety module connection

Neckpad : ECJ3900

- **Function**

A class I medical device that when connected to an orthodontic appliance, serves as a support for pulling.
Max 2,5kg safety module release force load capacity per side

- **Indications**

To be used during orthodontic treatments, based on the doctor's recommendation

- **Contraindications**

Allergic reactions while wearing the headcap during treatments.

Material : 100% cotton


Wash at 40°C with a general detergent used for cotton.

Discoloration risk : never wash together with other color headcap.

The product must not be reused if it is damaged. Any serious incident involving the device should be reported to the manufacturer and to the competent authority of the country in which the user resides.

Usage

safety module connection