



ENERGY CHAIN™

Endurance To Go The Distance


























Energy Chain™ outperforms competitors start to finish. Its unique patented thermosetting resin formula provides superior elastic memory following initial stretch over competitor's plastic chain.

University and clinical research demonstrate that RMO®'s Energy Chain delivers sustained light continuous force over a longer period than other available elastic chains. The tough patented material is stain-resistant, tasteless, non-toxic, and latex free.

Features and Benefits:

- Package contains a spool that is 7'6" (2.3m), Clear and Gray available in 15'42"(4.7m)
- Patented formula provides superior shape memory following initial stretch
- Close and open space faster with consistent force up to eight weeks
- Reduce appointment intervals and save valuable chair time
- Less stress decay and less elongation over time compared to other available elastic chains
- Spacing variation between modules allows proper selection





Rings size	Small	Small	Small	Small	Large	
Space between rings	 Closed	 Adjoined	 Reduced	 Medium	 Large	
Space between centers	2,8 mm	3 mm	3,8 mm	4 mm	5,1 mm	
Spool length	2,30 m	2,30 m	2,30 m	2,30 m	2,30 m	
Color						
Clear		J00120	J00202	J00269	J00121	J00122
Gray		J00166	J00203	J00270	J00167	J00168
Blue		J00223	J00205	J00271	J00224	J00281
Green		J00227	J00206	J00272	J00228	J00282
Red		J00231	J00204	J00273	J00232	J00280
Purple		J00235	J00207	J00274	J00236	J00283
Teal		J00136	J00214	J00137	J00138	J00290
Pearl Blue		J00139	J00215	J00140	J00141	J00291
Orange		J00239	J00208	J00275	J00240	J00284
Porsche Red		J00261	J00209	J00276	J00262	J00285
Black		J00263	J00210	J00277	J00264	J00286
Light Pink		J00265	J00211	J00278	J00266	J00287
Yellow		J00267	J00212	J00279	J00268	J00288
Ivory		J00142	J00213	J00143	J00144	J00289
Pink		J00145	J00216	J00146	J00147	J00292

ENERGY CHAIN™

Endurance To Go The Distance



Spool length	4,70 m	4,70 m	4,70 m	4,70 m	4,70 m	
Color						
Clear		J00293	J00295	J00297	J00299	J00301
Gray		J00294	J00296	J00298	J00300	J00302

Consider other products made of the same quality material as Energy Chain™:

Energy Thread™ Elastomeric Thread

RMO® Energy Thread™ is produced using the same material found in Energy Chain™. Energy Thread™ is extruded in two cross-section square sizes - .018" x .018" and .026" x .026." By producing Energy Thread™ in a square configuration, it ties snugly and eliminates knot slippage common with round thread.

- 20% more mass than round thread for enhanced force delivery
- Resists swelling, staining and discoloration
- Smooth surface and hydrophobic properties prevent plaque build-up better than round textured thread
- 25-feet per spool in clear or gray

Size	Color	Package Contains	Order Number
.025" (0.635mm), 7,5m	Clear	1 spool	J00386
.025" (0.635mm), 7,5m	Gray	1 spool	J00387
.030" (0.762mm), 7,5m	Clear	1 spool	J00388
.030" (0.762mm), 7,5m	Gray	1 spool	J00389
.018" x .018" (0.457mm x 0.457mm), 7,5m	Clear	1 spool	J00193
.018" x .018" (0.457mm x 0.457mm), 7,5m	Gray	1 spool	J00195
.026" x .026" (0.660mm x 0.660mm), 7,5m	Clear	1 spool	J00194
.026" x .026" (0.660mm x 0.660mm), 7,5m	Gray	1 spool	J00196



Tough-O® Light Energy Ligatures

Just like the original Tough-O®s, these ligatures maintain a higher force level between once visits compared to standard ligatures, delivering treatment efficiency by improving slot prescription expression and decreasing arch wire slop.

- Same great material as Energy Chain™
- Formulated for easy ligation
- Decreased elasticity deformation
- High force retention

Package Contains	Order Number
30 sticks (20ligs/stick)	J00159

