





















## ENERGY CHAIN™ F.M. Ringlet

Outperforms all other elastic chains

### NEW SIZE 2020

- Patented formula provides light continuous force for weeks
- Independently tested and clinically proven performance can reduce appointment intervals and save valuable chair time
- Less stress decay and less elongation over time compared to other elastic chains available
- Stain resistant and latex-free
- All Energy Chain™ colors perform similarly to Gray and Clear

|                       |  | NEW  |   |  |  |         |
|-----------------------|--|--|---|--|--|---------|
|                       | Closed  | Adjoined  | Reduced  | Medium  | Large  |         |
| Rings size            | Small  | Small  | Small   | Small  | Large  |         |
| Space between rings   | Closed   | Adjoined   | 1 mm  | 1,2 mm   | 1,2 mm   |         |
| Space between centers | 2,8 mm   | 3 mm   | 3,8 mm  | 4 mm   | 5,1 mm   |         |
| Spool length          | 4 m  | 4 m  | 4,70 m  | 4,70 m   | 4,70 m   |         |
| Couleur               |  |  |   |  |  |         |
| Clear                 |         | ECM0693  | <b>ECM0600</b>  | ECM0650  | ECM0691  | ECM0692 |
| Gray                  |        | ECM0698  | <b>ECM0610</b>  | ECM0660  | ECM0696  | ECM0697 |
| Red                   |       | ECM0661  | <b>ECM0601</b>  | ECM0651  | ECM0671  | ECM0681 |
| Blue                  |       | ECM0662  | <b>ECM0602</b>  | ECM0652  | ECM0672  | ECM0682 |
| Green                 |       | ECM0663  | <b>ECM0603</b>  | ECM0653  | ECM0673  | ECM0683 |
| Violet                |       | ECM0664  | <b>ECM0604</b>  | ECM0654  | ECM0674  | ECM0684 |
| Orange                |       | ECM0665  | <b>ECM0605</b>  | ECM0655  | ECM0675  | ECM0685 |
| Porsche Red           |       | ECM0666  | <b>ECM0606</b>  | ECM0656  | ECM0676  | ECM0686 |
| Black                 |       | ECM0667  | <b>ECM0607</b>  | ECM0657  | ECM0677  | ECM0687 |
| Light Pink            |       | ECM0668  | <b>ECM0608</b>  | ECM0658  | ECM0678  | ECM0688 |
| Yellow                |       | ECM0669  | <b>ECM0609</b>  | ECM0659  | ECM0679  | ECM0689 |
| Ivory                 |       | ECM0620  | <b>ECM0611</b>  | ECM0621  | ECM0622  | ECM0623 |
| Teal                  |       | ECM0624  | <b>ECM0612</b>  | ECM0625  | ECM0626  | ECM0627 |
| Pearl Blue            |       | ECM0628  | <b>ECM0613</b>  | ECM0629  | ECM0630  | ECM0631 |
| Pink                  |       | ECM0632  | <b>ECM0614</b>  | ECM0633  | ECM0634  | ECM0635 |

