


Altitude SL® M Brackets

Emphasizing more on inter-bracket distance, lower frictional resistance, and smoothness. We have created an extra small profile, rounded slot edge, and a fully rounded contour design.

Altitude SL® M, the self-ligating bracket that gives you more control and comfort from start to finish.

MATERIAL: 17-4 PH Stainless Steel
BASE: 304 Stainless Steel 80-gauge Mesh
PRESCRIPTION: Roth
SLOT SIZE: .018 & .022



FEATURES	BENEFITS
MIM technology	High level precision
Fully rounded contour design	Extra comfort
Extra small profile design	Increased inter-bracket distance
Patented locking mechanism	Opening with a regular probe and easy closing with a finger
Sheet metal slot and rounded slot edge	Reduced frictional resistance
80-gauge mesh base	Superior bonding strength
Laser mark on the base	Easy bracket identification

Roth Prescription - Maxillary Brackets

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central	12°	5°	0°		R	M14300	M14400
	12°	5°	0°		L	M14301	M14401
Lateral	8°	9°	0°		R	M14302	M14402
	8°	9°	0°		L	M14303	M14403
Cuspid	-2°	10°	0°	Distal	R	M14306	M14406
	-2°	10°	0°	Distal	L	M14307	M14407
1st/2nd Bicuspid	-7°	0°	0°	Distal	R	M14310	M14410
	-7°	0°	0°	Distal	L	M14311	M14411



Altitude SL M_29/08/2019

Prescription Roth - Brackets Mandibulaires

Tooth	Torque	Angulation	Offset	Hook	R/L	.018	.022
Central/ Lateral	-1°	0°	0°		R/L	M14320	M14420
Cuspid	-11°	7°	0°	Distal	R	M14326	M14426
	-11°	7°	0°	Distal	L	M14327	M14427
1st Bicuspid	-17°	0°	0°	Distal	R	M14330	M14430
	-17°	0°	0°	Distal	L	M14331	M14431
2nd Bicuspid	-22°	0°	0°	Distal	R	M14334	M14434
	-22°	0°	0°	Distal	L	M14335	M14435

Roth Prescription - Kit

Description	.018	.022
Kit with hook on cuspid/bicuspid, 5-5(20 brackets)	KM14302	KM14402